

# **An example of packaging deep learning model "AIGenMols" for MolAICal**

Qifeng Bai, Shuoyan Tan  
Email: molaical@yeah.net  
Homepage: <https://molaical.github.io>  
Lanzhou University  
Lanzhou, Gansu 730000, P. R. China

## 1. Introduction

Sometimes, the "AIGenMols" of MolAICal (<https://doi.org/10.1093/bib/bbaa161>) have not good compatibility in Linux operating system. In this case, it needs to generate the binary "AIGenMols" again. In this tutorial, the ORGAN (<https://github.com/gablg1/ORGAN>) is chosen for the installation of "AIGenMols". Besides, you can build the "AIGenMols" from your trained deep learning model.

## 2. Materials

### 2.1. Software requirement

- 1) MolAICal: <https://molaical.github.io>
- 2) Anaconda: <https://www.anaconda.com>

You can choose the free version of individual edition of Anaconda. The version of Anaconda should be chosen based on Python 3.x rather than Python 2.x, etc.

### 2.2. Example files

- 1) All the necessary tutorial files are downloaded from:  
<https://github.com/MolAICal/specialtopic/tree/master/012-AIGenMols>

## 3. Procedure

- 3.1. Install Anconda (<https://www.anaconda.com/distribution>).

Make sure installation completely and environment variables are set rightly.

- 3.2. Download "spec-list.txt" and create an independent environment

```
#> conda create --name AIGen --file spec-list.txt
```

```
#> conda activate AIGen
```

**Note:** If some errors happen like: *CondaHTTPError: HTTP 404 NOT FOUND for url <<https://conda.anaconda.org/conda-forge/linux-64/xxxx-xxxx>>*, you can delete corresponding soft package in "spec-list.txt". Then repeat above steps. The missing soft package can be installed without version number via similar below step:

```
#> conda install -c conda-forge xxxx
```

- 3.3. Install the necessary libraries

```
#> conda install -c conda-forge openblas
```

```
#> conda install -c conda-forge openssl
```

```
#> pip install editdistance==0.3.1
```

- 3.4. Finally, produce binary "deep learning model".

Decompress "organ.zip"

```
#> unzip organ.zip
```

```
#> cd organ
```

```
#>          pyinstaller          --add-data="organ/NP_score.pkl.gz:organ"          --add-
data="organ/SA_score.pkl.gz:organ"          --add-data="organ/data/FDA-H.csv:organ/data"          --add-
data="organ/checkpoints/FDA-H/FDA-H_99.ckpt.data-00000-of-00001:organ/checkpoints/FDA-
H"          --add-data="organ/checkpoints/FDA-H/FDA-H_99.ckpt.index:organ/checkpoints/FDA-H" --
add-data="organ/checkpoints/FDA-H/FDA-H_99.ckpt.meta:organ/checkpoints/FDA-H"          --add-
data="organ/data/FDA1884.csv:organ/data"          --add-
data="organ/checkpoints/FDA1884/FDA1884_119.ckpt.data-00000-of-
00001:organ/checkpoints/FDA1884"          --add-
data="organ/checkpoints/FDA1884/FDA1884_119.ckpt.index:organ/checkpoints/FDA1884" --
add-data="organ/checkpoints/FDA1884/FDA1884_119.ckpt.meta:organ/checkpoints/FDA1884"
--add-data="organ/data/zinc.csv:organ/data"          --add-
data="organ/checkpoints/ZINC/ZINC_99.ckpt.data-00000-of-00001:organ/checkpoints/ZINC" -
-add-data="organ/checkpoints/ZINC/ZINC_99.ckpt.index:organ/checkpoints/ZINC"          --add-
data="organ/checkpoints/ZINC/ZINC_99.ckpt.meta:organ/checkpoints/ZINC" -F main.py
```

Open folder “organ/dist”

```
#> cd dist
```

You will see file named “main”, rename it as “AIGenMols”

```
#> mv main AIGenMols
```

```
#> chmod +x AIGenMols
```

The deep binary deep learning model named “AIGenMols” is produced completely. Put it into MolAICal software folder directly.

### Analysis of install process:

```
-----
--add-data="organ/data/zinc.csv:organ/data"          --add-
data="organ/checkpoints/ZINC/ZINC_99.ckpt.data-00000-of-00001:organ/checkpoints/ZINC" -
-add-data="organ/checkpoints/ZINC/ZINC_99.ckpt.index:organ/checkpoints/ZINC"          --add-
data="organ/checkpoints/ZINC/ZINC_99.ckpt.meta:organ/checkpoints/ZINC" -F main.py
-----
```

You can replace above data by trained deep learning model.

### Notice:

1). When you install Anaconda, the install path of Anaconda should have no blank space for packaging “AIGenMols”. For example, /home/test install/anaconda. The blank space between the characters “test” and “install” may lead to the failure of the packaging process.

2). The command usage of “AIGenMols” should be like as below:

```
AIGenMols -f test.txt -s FDAFrag -n 100
```

```
AIGenMols -f test.txt -s ZINCMol -n 100
```

- f:** represents file name which contains generated results of molecules with SMILES format
- s:** represents selection of deep learning model
- n:** represents the number of generated molecules

This is just an example for packaging “AIGenMols”. You can use the similar way for your trained model. We will update deep learning model.